





SOP: TaifaCare Mortuary Management Module – Compartment Setup User Guide


(Last updated: August 2025)

Tasks:	Mortuary and Compartment Allocation Management			
Objective:	To guide KenyaEMR users through the process of setting up and managing mortuary compartment using the Bed Management module.			
Who:	EMR Admin User/HRIO/Mortician/Pathologist			
Required Materials:	Username, password, computer installed with TaifaCare version 19.2.0 and above			
Step	Action			
Login to KenyaEMR 3.X On successfully accessing the TaifaCare system you should be navigated to the login page, Authenticate entry by providing username and password on the fields, then click "Login" button.	TaifaCære Username Powered by A distribution of KenyaEMR CopenMRS CopenMRS CopenMRS			
















Enter the name of Tag Location the Mortuary you would like to create Mortuary Location Visit Location Login Location Admission Location e.g. **'KNH** Mortuary', Ensure 4 × Select a tag you add the correct Admission Location location Tag i.e. ✓ Login Location Admission ✓ Visit Location Location, Amenity Ward **Mortuary Location** and Visit Location then save the form. Cancel
















type e.g. **Std Comp** to denote when adding a Compartment that it will be a basic compartment for accommodating a deceased patient) then click on **Save** to submit the form









3. You should get a success pop up message and see the bed type added.
Please note you can also **Edit** the bed type details in case of any changes you would like to make in the future by clicking on the 'Pencil'

looking icon under

Actions



Bed Type


© Kayole II Sub-District Hospital · ☐ Today, 08:31

Bed Type			Add Bed Type +
Name	Display Name	Description	Actions
Standard	Standard	Inpatient	<u>@</u>
Deluxe	Deluxe		<u>@</u>
Cot	Cot	Cot for children	<u>@</u>
Cot	Cot	Cot for Children	<u>@</u>
ken	delux	kd	<u>@</u>
01-GYNAE	GYNAE		<u>@</u>
Cots	Paeds	Cots	<u>@</u>
ОТ	ОТ	ORTHOPEDIC	<u>@</u>
BBB	Paeds	R	<u>@</u>
surgy12	surg 2	standard	<u>@</u>
Std COMP	Std Comp	Standard Compartment For the deceased body	<u>@</u>
Items per page: 10 ∨ 1–10 of 11 items			
















You can also edit TaifaCære | ® Kayole II Sub-District Hospital Bed Management the bed details by clicking on the Ward Allocation Ward Allocation 'Pencil' icon under Edit bed Actions then make changes before Ward Alloc COMP001 saving. Bed ID Standard Compartment For the deceased body Bed row MW-003 2 Bed column 5 Location surg 1 Morgue × V Occupied Status Available Bed type Std COMP MAT-001









THE END