





### **How to Record Vitals in TaifaCare**

Tasks:	To Record vitals on TaifaCare	
Objective:	Submit Taifa Care Triage Form	
Who:	EMR user/Triage Nurse/Clinician	
Required Materials:	Username, password, computer installed with TaifaCare latest version	
Version:	Taifa Care Version 19.2.0 and above	
Last Updated:	Aug 2025	

Step	Action	Screen shots
Introduction	This document provides a detailed workflow on how to Record vitals in the Taifa Care system. The guide is accompanied by relevant screenshot images for illustration purposes.	








Log into Taifa Care	The user is required to gain access into the Taifa Care system by supplying correct username and password.  - Enter Username  - Enter Password  - Click "Login"	TaifaCæe  Username  Password  Log in  A distribution of
KenyaEMR landing page	After a successful login, the Taifa Care landing page (Active visits) will be opened. Click Home to navigate to the TaifaCare Homepage.	(♣) KenyaEMR ○ OpenMRS ○ OpenMRS
















## Open the Triage form and Recording vitals

After checking-in the client, proceed to click on the 'Record vitals' linkage on the patient profile (as highlighted) and open the Triage form.

Or after checking-in the client Navigate to the right panel, click on the clinical forms to open the forms, then search for the triage form on the search field.


Click on the Triage form to record the vitals.



























#### **Recording Anthropometrics**

#### 3. Anthropometrics and other Recordings

- Record Weight(kg) and Height(cm) in Numeric
- N/B; You can use the previous value for Height.
- BMI (kg/m2) is auto calculated.
- Nutritional status is being auto populated.

3.Do pregnancy screening if eligible.

**NB:** This is only availed for women.










# **Completing Vital recording**

- 4.Recording Triage Notes
  - Record triage notes if any in the text field provided on the 'Triage Notes'
  - BMI (kg/m2) is auto calculated.
  - Nutritional status is being auto populated.


5.Click **'Save and close'** to save the form





























THE END